





INSTRUCTION OF USE: Seperator CC / Seperator HC

Material information and manufacturer data

Name, trade name	Seperator CC / Seperator HC				<input type="checkbox"/> Medical Device
Classification	Classification:	-	Labelling:		None
Data	Article No.:	F73, F74		Basis UDI:	763029670073P, 763029670074P
Variants	Powder	Liquid		Colours	
	-	100ml, 250ml, 500ml, 5.000ml		-	
Allowed equipment	-				
Manufacturer data	SRN:		LI-MF-000000633		
	Anschrift:		NOVODENT Ets., Industriestrasse 32, 9487 Bendorf, Liechtenstein		
General information	Contamination of the insulating fluid by gypsum residues, metal-clad brushes, or similar materials should be avoided.				
Intended purpose	Insulation of gypsum models in plastics engineering.				
Indication	Suitable for insulating gypsum models in plastics engineering. The insulating film ensures smooth surfaces on the plastic and ensures flawless separation between plastic and gypsum as well as between gypsum and gypsum.				
Contraindication	Not suitable for use in the mouth. Not approved as a medical device.				

Use of the product

Properties
The separator forms a very thin, smooth, precise, and uniform insulation. It ensures flawless separation of plastic and gypsum. The separator can be applied to cold or warm gypsum.

Environment of use	Shelf life and storage instructions
Applicable at room temperature.	Avoid direct sunlight. Always keep containers closed. The shelf life is 2 Years. Do not use the product after the expiry date. Keep out of reach of children. Store at room temperature min. + 2°C, max. 25°C.
Material composition	
Liquid: Alginate hydrogel	
Skin protection	Alginate hydrogel is skin-friendly. Individuals with allergic sensitivities to the ingredients may wear latex gloves for protection.
Transport	For regulations on the transport of the raw material, please refer to the Seperator CC / Seperator HC safety data sheets.
Specifications	No specific product specifications available. Seperator CC / Seperator HC
















Regulatory information	-
Disposal according to national regulations.	

Instructions for use

Processing	<p>Preparation and use Before isolating, all dried gypsum models must be thoroughly moistened (5-10 minutes). This prevents air from getting trapped under pressure in the soft plastic dough during polymerization.</p> <p>After the wax has been boiled out, the required amount of separator is placed in a container and applied 1- to 2-times to the surface to be insulated with a clean brush. The separator can be smoothed or dried with a gentle stream of compressed air. It is important to ensure that no insulation fluid remains in grooves and depressions and that no plastic teeth or clasps come into contact with the insulation. Once the insulating film is completely dry, the plastic dough can be poured immediately.</p>
-------------------	---

<p>Liability notice</p> <p>The information is based on the current state of our knowledge, but does not constitute a guarantee of product properties and does not establish a contractual legal relationship. The product is normally delivered stabilised. However, it may polymerise under heat development if the storage time and/or storage temperature are significantly exceeded. No liability is accepted for damage resulting from failure to observe the instructions for use and the area of application. This also applies if the product is mixed or processed together with products from competitors.</p>

Contact details of the manufacturer			
Address:	NOVODENT Ets., Industriestrasse 32, 9487 Bendern, Liechtenstein		
Contact:	Phonel.: +423 377 1500	Fax: +423 377 1509	E-Mail: frontoffice@novodent.com
Websites::	Webshop: www.novodent.dental	Main page:	www.novodent.com

Note about the symbols ⓘ			
	<input checked="" type="checkbox"/> Lot/Batch No.		<input checked="" type="checkbox"/> Manufacturer
	<input checked="" type="checkbox"/> Expiry date		<input checked="" type="checkbox"/> Storage temperatures
	<input checked="" type="checkbox"/> Articlenuumber		<input type="checkbox"/> Protection from sunlight
	<input type="checkbox"/> Attention (GHS)		<input type="checkbox"/> Danger Flammable (GHS)
	<input checked="" type="checkbox"/> Single Registration Number (EUDAMED)		<input checked="" type="checkbox"/> Unique Device Identification 
			<input checked="" type="checkbox"/> Observe operating instructions
			<input type="checkbox"/> Medical Device *
			<input checked="" type="checkbox"/> Do not reuse
			<input checked="" type="checkbox"/> MP- European conformity *

Symbols used on the labelling of this product are marked with . *Not applicable for non-medical products